



## Optos Suzuka DICOM Conformance Statement


Company Name: **OPTOS**

Product Name: **Suzuka**

Version: 017

Document Number: **D108878**

Date: **20/12/2021**

<p>Author</p>	 <p>Colin Keillo Dec 20 2021 11:47 AM</p>
<p>Product Management</p>	 <p>Iain Gourlay Dec 20 2021 12:00 PM</p>
<p>Quality Assurance</p>	 <p>Matthew Martin Dec 20 2021 11:58 AM</p>

This document, and the information contained herein, is the property of Optos PLC. It is supplied commercially in confidence without liability for errors or omissions and must not be reproduced, used or disclosed, in part or in whole, without prior written permission.

© Copyright 2021 All Rights Reserved, Optos PLC.

# 1 Conformance Statement Overview

This document describes the DICOM conformance of Optos devices, P200T, P200DTx, P200TE and P200TxE when used with Optos Suzuka image processing software. The system implements the necessary DICOM services to save acquired ophthalmic images to a network storage device. The system conforms to the DICOM standard to allow the sharing of medical information with other digital imaging systems.

Table 1 provides an overview of the network services that are supported by the OPTOS devices with OPTOS Suzuka Image Processing Software.

**Table 1: Network Services**

SOP Class	User Of Service (SCU)	Provider of Service (SCP)
<b>Transfer</b>		
Grayscale Softcopy Presentation State Storage	Yes	No
Ophthalmic Photography 8-bit Image Storage	Yes	No
Ophthalmic Tomography Image Storage	Yes	No
Widefield Ophthalmic Photography Stereographic Projection Image Storage	Yes	No
Stereometric Relationship Storage	Yes	No

Table 2 provides an overview of the Media Storage Application Profiles supported by the Software.

**Table 2: Media Services**

Media Storage Application Profile	Write Files (FSC or FSU)	Read Files (FSR)
*	No	No

## 2 Table of Contents

<b>1</b>	<b>Conformance Statement Overview</b> .....	<b>2</b>
<b>2</b>	<b>Table of Contents</b> .....	<b>3</b>
<b>3</b>	<b>Introduction</b> .....	<b>5</b>
3.1	Revision History .....	5
3.2	Audience .....	6
3.3	Remarks .....	7
3.4	Terms and Definitions .....	8
3.5	Basics of DICOM Communication .....	10
3.6	Abbreviations .....	10
3.7	References .....	11
<b>4</b>	<b>Networking</b> .....	<b>12</b>
4.1	Implementation Model .....	12
4.1.1	Application Data Flow Diagram .....	12
4.1.2	Functional Definition of Storage Application Entity .....	12
4.1.3	Sequencing of Real-World Activities .....	13
4.2	Implementation Model – Modality Worklist .....	13
4.2.1	Application Data Flow Diagram .....	13
4.2.2	Functional Definitions of Modality Work List Application Entity .....	13
4.2.3	Sequencing of Real-World Activities .....	13
4.3	Implementation Model - Storage .....	14
4.3.1	Functional Definitions of Storage Application Entities .....	14
4.3.2	Sequencing of Real World Activities .....	14
4.4	Application Entity Specifications .....	14
4.4.1	Storage Application Entity Specification .....	14
4.4.1.1	SOP Classes .....	14
4.4.1.2	Association Policies .....	15
4.4.1.3	Association Initiation Policies .....	15
4.4.1.4	Association Acceptance Policy .....	19
4.5	Network interfaces .....	25
4.5.1	Physical Network Interface .....	25
4.5.2	Additional Protocols .....	25
4.5.3	IPv4 and IPv6 Support .....	25
4.6	Configuration .....	26
4.6.1.1	Application Entity Title/Presentation Address Mapping .....	26
4.6.1.2	Remote Application Entity Title / Presentation Address Mapping .....	26
4.6.1.3	Parameters .....	26
<b>5</b>	<b>Support for Extended Character Sets</b> .....	<b>28</b>
5.1	Character sets in use .....	28
5.1.1	Character set configuration .....	28
<b>6</b>	<b>Security</b> .....	<b>28</b>
6.1	Conformance .....	28
<b>7</b>	<b>Annexes</b> .....	<b>29</b>
7.1	IOD Contents .....	29
7.1.1	Created SOP Instances .....	29
7.1.1.1	Ophthalmic Photography 8-Bit Image IOD .....	30
7.1.1.2	Widefield Ophthalmic Photography Stereographic Projection IOD .....	30
7.1.1.3	Ophthalmic Tomography Image IOD .....	31
7.1.1.4	Grayscale Softcopy Presentation State IOD .....	34
7.1.1.5	Stereometric Relationship IOD .....	35
7.1.1.6	Common Modules .....	36

---

7.1.1.7	Ophthalmic Photography 8 Bit Image Modules .....	40
7.1.1.8	Widefield Ophthalmic Photography Stereographic Projection Image Modules ....	46
7.1.1.9	Ophthalmic Tomography Image Modules.....	53
7.1.1.10	Grayscale Softcopy Presentation State Modules .....	63
7.1.1.11	Stereometric Relationship Modules .....	66
7.2	Usage of Attributes from received IOD's .....	67
7.3	Attribute Mapping.....	67
7.4	Coerced / Modified Fields .....	67
7.5	Data Dictionary of Private Attributes .....	67
7.6	Coded terminology and Templates .....	70
7.7	Grayscale Image Consistency .....	70
7.8	Standard/Extended/Specialized/Private SOP Classes .....	70
7.9	Private Transfer Syntaxes .....	70

## 3 Introduction

### 3.1 Revision History

Table 3: Revision History

Revision	Date	Author	Notes
001	19 <sup>th</sup> April 2015	Anilasha Varkey	First edition- DCS for P200TE
002	1 <sup>st</sup> December 2015	Gavin Falconer	Second edition DCS for P200TE
002	15 <sup>th</sup> Feb 2016	Gavin Falconer	Corrections to Second edition DCS for P200TE: Table 54 and table 58 for OPT SOP instances; Tables 59, 60 and 62 for GSPS SOP instances.
003	17 <sup>th</sup> May 2016	Gavin Falconer	Third edition DCS for P200TE, with Suzuka release v1.0.1. Added Private Ophthalmic Frame Location Functional Group Macro for OPT SOP Instances. Corrections to VR for several OPT SOP instance tags.
004	29 <sup>th</sup> Sept 2016	Gavin Falconer	Fourth edition DCS for P200TE with Suzuka release v1.1.0.
005	8 <sup>th</sup> November 2016	Gavin Falconer	Added private tags for ProView sharpening.
005	13 <sup>th</sup> December 2016	Anilasha Varkey	Added private tags for OCT scan metrics
006	16 <sup>th</sup> February 2017	Anilasha Varkey	Bug 8601 fix- corrected the 'Ophthalmic Image Orientation' tag value for RNFL image Bug 8606 fix- correct dicom tag number for ICC Profile.
007	27 <sup>th</sup> February 2017	James Howat	Bug 8618 fix – correct dicom VRs for scan outcome labels.
008	7 <sup>th</sup> September 2017	Steven Reid	Bug 8803 fix – correct scan type names. Align private tag definitions with Vantage.
009	12 <sup>th</sup> April 2018	James Howat Steven Reid	Bug 9252 fix – align General Equipment Manufacturer tag value with Vantage. Added private tag for OCT fovea location and OCT segmentation outcomes. Added private tag usages for General Projection.

010	7 <sup>th</sup> November 2018	James Howat	Updated private tags for OCT fovea location, disc location and 3D cup and disc perimeter. Updated OPT Equipment values for P200TE.
011	25 <sup>th</sup> January 2019	Michael Irvine James Howat	Added private tag for Optos Procedure Description Merged content for P200TxE Added Stereometric Relationship Storage SOP Class
012	20 <sup>th</sup> March 2019	Michael Irvine	Added entries to Table 52 to support linking of ONH and derived RNFL images
013	11 <sup>th</sup> May 2020	Gary Young	Added entries into Table 52 to support linking of OCT scans to predecessors Added entries into tables 61 and 71 to support follow up scans Added details of follow up data that is now being stored in the dicoms in section 4.2.1.3.1.1 Description and sequencing of activities
014	17 <sup>th</sup> July 2020	Mike Irvine Steven Reid	Added Patient Comments entry into Table 18 Update unprojected SLO image private tag
015	19 <sup>th</sup> April 2021	Steven Reid Colin Keillor Mike Irvine	Add Segmentation Fovea Detection Confidence and Outcome into Table 61 and Table 71 Updated Study Instance UID tags value and source details. Add Operator's Name entry into Table 20 Updated 4.2.1 and 4.4.1.2 to cover C-Store
016	07 <sup>th</sup> September 2021	Colin Keillor	Added sections 4.2 and 4.4.1.4.1 to 4.4.1.4.1.3 detailing MWL support.
017	20 <sup>th</sup> December 2021	Colin Keillor	Updated section 1 adding list of supported Optos devices. Correction of table numbers throughout document.

### 3.2 Audience

This document is written for the people that need to understand how the OPTOS Device and Suzuka Image Processing Software will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product.

This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

### 3.3 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between the OPTOS Software and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard. DICOM by itself does not guarantee interoperability.

The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality. This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

---

## 3.4 Terms and Definitions

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

**Abstract Syntax** – the information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.

**Application Entity (AE)** – an end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.

**Application Entity Title** – the externally known name of an Application Entity, used to identify a DICOM application to other DICOM applications on the network.

**Application Context** – the specification of the type of communication used between Application Entities. Example: DICOM network protocol.

**Association** – a network communication channel set up between Application Entities.

**Attribute** – a unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

**Information Object Definition (IOD)** – the specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as:

- Mandatory (Type 1)
- Required but possibly empty (Type 2)
- Optional (Type 3)

There may be also be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.

**Joint Photographic Experts Group (JPEG)** – a set of standardized image compression techniques, available for use by DICOM applications.

**Media Application Profile** – the specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs)

**Module** – a set of Attributes within an Information Object Definition that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.

**Negotiation** – first phase of Association establishment that allows Application Entities to agree on the types of data to be exchanged and how that data will be encoded.

**Presentation Context** – the set of DICOM network services used over an Association, as negotiated between Application Entities; includes Abstract Syntaxes and Transfer Syntaxes.

---

**Protocol Data Unit (PDU)** – a packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.

**Security Profile** – a set of mechanisms, such as encryption, user authentication, or digital signatures, used by an Application Entity to ensure confidentiality, integrity, and/or availability of exchanged DICOM data

**Service Class Provider (SCP)** – role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity (Service Class User). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).

**Service Class User (SCU)** – role of an Application Entity that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)

**Service/Object Pair (SOP) Class** – the specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.

**Service/Object Pair (SOP) Instance** – an information object; a specific occurrence of information exchanged in a SOP Class. Examples: a specific x-ray image.

**Tag** – a 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the “group” and the “element”. If the “group” number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

**Transfer Syntax** – the encoding used for exchange of DICOM information objects and messages. Examples: JPEG compressed (images), little-endian explicit value representation.

**Unique Identifier (UID)** – a globally unique “dotted decimal” string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.

**Value Representation (VR)** – the format type of an individual DICOM data element, such as text, an integer, a person’s name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

**C-ECHO** - DICOM service used to verify end-to-end communications with a peer DICOM-service-user

**C-STORE** - DICOM service for storing objects, e.g. images

**Modality Worklist** - DICOM service to inform the modality of scheduled examinations

### 3.5 Basics of DICOM Communication

This section describes terminology used in this Conformance Statement for the non-specialist.

The key terms used in the Conformance Statement are highlighted in italics below. This section is not a substitute for training about DICOM, and it makes many simplifications about the meanings of DICOM terms.

Two Application Entities (devices) that want to communicate with each other over a network using DICOM protocol must first agree on several things during an initial network “handshake”. One of the two devices must initiate an Association (a connection to the other device), and ask if specific services, information, and encoding can be supported by the other device (Negotiation).

DICOM specifies a number of network services and types of information objects, each of which is called an Abstract Syntax for the Negotiation. DICOM also specifies a variety of methods for encoding data, denoted Transfer Syntaxes. The Negotiation allows the initiating Application Entity to propose combinations of Abstract Syntax and Transfer Syntax to be used on the Association; these combinations are called Presentation Contexts. The receiving Application Entity accepts the Presentation Contexts it supports.

For each Presentation Context, the Association Negotiation also allows the devices to agree on Roles – which one is the Service Class User (SCU - client) and which is the Service Class Provider (SCP - server). Normally the device initiating the connection is the SCU, i.e., the client system calls the server, but not always.

The Association Negotiation finally enables exchange of maximum network packet (PDU) size, security information, and network service options (called Extended Negotiation information). The Application Entities, having negotiated the Association parameters, may now commence exchanging data. Common data exchanges include queries for worklists and lists of stored images, transfer of image objects and analyses (structured reports), and sending images to film printers. Each exchangeable unit of data is formatted by the sender in accordance with the appropriate Information Object Definition, and sent using the negotiated Transfer Syntax. There is a Default Transfer Syntax that all systems must accept, but it may not be the most efficient for some use cases. Each transfer is explicitly acknowledged by the receiver with a Response Status indicating success, failure, or that query or retrieve operations are still in process.

### 3.6 Abbreviations

AE	Application Entity
AET	Application Entity Title
CD-R	Compact Disk Recordable
DICOM	Digital Imaging and Communications in Medicine
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
IHE	Integrating the Healthcare Enterprise
IOD	Information Object Definition
IS	Information System (HIS, RIS, PACS within hospital)
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standards
JPEG	Joint Photographic Experts Group
LUT	Look-up Table
MWL	Modality Worklist

OP	Ophthalmic Photography
OSI	Open Systems Interconnection
PACS	Picture Archiving and Communication System
PDU	Protocol Data Unit
SC	Secondary Capture
SCP	Service Class Provider
SCU	Service Class User
SOP	Service-Object Pair
TCP/IP	Transmission Control Protocol/Internet Protocol
UID	Unique Identifier
VR	Value Representation
WSD	Workstation

### 3.7 References

NEMA PS3 Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>

## 4 Networking


The supported OPTOS devices with OPTOS Suzuka Image Processing Software allows for the following DICOM functionality:

- Exporting DICOM images to a DICOM SCP

### 4.1 Implementation Model

#### 4.1.1 Application Data Flow Diagram

This DICOM Conformance statement specifies the behavior and functionality of the OPTOS Suzuka Image Processing Software.



The Storage Application Entity sends images, Presentation States and, where applicable, Stereometric Relationship objects, to a remote AE. It is associated with the local real-world activity "Save exam". "Save exam" is performed automatically when a scan head operator completes a patient exam.

The images associated with the real-world activity "Save exam" may include one or more:

- Widefield Ophthalmic Photography images from a Scanning Laser Ophthalmoscope (SLO). A Stereographic Projection transform will have been applied to these images.
- Optical Coherence Tomography images, where each image comprises a collection of OCT image frames, accompanied by an SLO index image. (An SLO index image will not have been subject to Stereographic Projection).

#### 4.1.2 Functional Definition of Storage Application Entity

The Storage Application Entity encodes image content acquired by an imaging device into appropriate DICOM SOP instances. The application then uses the DICOM protocol to transfer the data across the network to the remote SCP.


### 4.1.3 Sequencing of Real-World Activities

Not applicable

## 4.2 Implementation Model – Modality Worklist

### 4.2.1 Application Data Flow Diagram

The worklist is queried by the AE and supplies the SCU (Optos device and Suzuka software) with the scheduled tasks that have to be performed on the modality. Optos device and Suzuka software supports the worklist service in the role of a SCU.



### 4.2.2 Functional Definitions of Modality Work List Application Entity

The modality worklist SCU requests the worklist SCP to perform a match to the keys specified in the C-FIND DIMSE service. The worklist SCP responds to the C-FIND query and sends scheduled imaging service requests and patient demographic information from the information system to the modality. Suzuka software then sends the result from the C-FIND service to the Optos device to initiate the examination.

### 4.2.3 Sequencing of Real-World Activities

In order to perform a worklist query, the user must search for the scheduled procedure using the search criteria (patient based query or broad query) on the Optos device. On receiving a call to the patient search API endpoint, Suzuka software sends a C-FIND query with the search keys to the configured SCP. The matching results are then received by the Optos Scanhead Device via the Suzuka software. The displayed results are based on the search keys entered in the query fields. After the task is completed, the user can perform a new query.

## 4.3 Implementation Model - Storage

### 4.3.1 Functional Definitions of Storage Application Entities

Assuming the system is correctly configured Suzuka worker application will create an association with the designated Storage SCP and transfer the DICOM image files. This is performed as a transparent process from the user's perspective; there is no user feedback to indicate that the transfer has occurred or any other status regarding the transfer. If a problem arises then the detail of the error is logged.

### 4.3.2 Sequencing of Real World Activities

The Suzuka software implements only the Storage AE, so real-world activities are limited to the storing of images. There is, therefore, no scope for sequencing or interaction of separate AE's.

## 4.4 Application Entity Specifications

### 4.4.1 Storage Application Entity Specification

The software can be configured to store DICOM images to a PACS / storage system via DIMSE C-STORE.

#### 4.4.1.1 SOP Classes

This Application Entity provides Standard or Standard Extended Conformance to the following DICOM SOP Classes:

**Table 4: SOP Classes for AE Storage**

SOP Class Name	SOP Class UID	SCU	SCP
Ophthalmic Photography 8 Bit Image Storage SOP Class	1.2.840.10008.5.1.4.1.1.77.1.5.1	Yes	No
Widefield Ophthalmic Photography Stereographic Projection Image Storage SOP Class	1.2.840.10008.5.1.4.1.1.77.1.5.5	Yes	No
Ophthalmic Tomography Image Storage SOP Class	1.2.840.10008.5.1.4.1.1.77.1.5.4	Yes	No
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1	Yes	No
Stereometric Relationship Storage SOP Class	1.2.840.10008.5.1.4.1.1.77.1.5.3	Yes	No

#### 4.4.1.2 Association Policies

##### 4.4.1.2.1 General

The DICOM standard application context name for DICOM 3.0 is used:

**Table 5: DICOM application context for AE Storage**

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

The software requests a maximum PDU size of 16 Kbytes (16384 bytes).

The software does not support extended negotiations.

The only supported network protocol is TCP/IP. The underlying Operating System provides the TCP/IP stack.

The Application Entity Title, the port number, and the host name or net address get defined in the Image Processing service configuration file.

##### 4.4.1.2.2 Number of Associations

The AE will attempt only one association establishment at a time.

**Table 6: Number of Associations Initiated for AE Storage**

Maximum number of simultaneous Associations	1
---------------------------------------------	---

##### 4.4.1.2.3 Asynchronous Nature

The AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

**Table 7: Asynchronous Nature as an SCU for AE Storage**

Maximum number of outstanding asynchronous transactions	1
---------------------------------------------------------	---

##### 4.4.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

**Table 8: DICOM Implementation Class and Version for AE Storage**

Implementation Class UID	1.2.276.0.7230010.3.0.3.6.2
Implementation Version Name	OFFIS_DCMTK_362

##### 4.4.1.3 Association Initiation Policies

The software will initiate an association as an outcome of the following real-world activity:

- Save exam initiated by the operator for acquired images

##### 4.4.1.3.1 Activity – Save exam

###### 4.4.1.3.1.1 Description and sequencing of activities

When a scan head operator completes a patient exam, the Storage Application Entity is automatically invoked to generate appropriate SOP instances and transfer the acquired images to a DICOM storage SCP.

Where an exam includes Widefield Ophthalmic Photography images, each image will be encoded as either a Widefield Ophthalmic Photography Stereographic Projection Image Storage SOP instance, or as an Ophthalmic Photography 8-Bit Image Storage SOP instance, depending on the value of the "Support Widefield Ophthalmic Photography Stereographic Projection SOP Class" configuration setting for the application.

In addition, where an exam includes Widefield Ophthalmic Photography images, a Grayscale Softcopy Presentation State SOP instance will be generated for the exam, containing display settings for each of the Widefield Ophthalmic Photography images.

Where an exam includes Optical Coherence Tomography images, each such image will result in an Ophthalmic Photography 8-bit Image Storage SOP instance for the SLO index image, and an Optical Coherence Tomography Image Storage SOP instance for the OCT frames.

Where an exam includes any pairs of either Widefield Ophthalmic Photography Stereographic Projection Image Storage SOP or Ophthalmic Photography 8-Bit Image Storage SOP instances, identified by the scan head as being two parts of a stereometric pair, a Stereometric Relationship Storage SOP instance will be generated, recording the relationship between all pairs.

When an exam includes information of a previous exam, then the temporal predecessor value will be recorded.

When an exam provides scan information relating to a previous exam then the information will be recorded in the OCT Module.

#### **4.4.1.3.1.2 Presentation Contexts**

The Storage Application Entity exports files within which the Presentation Contexts are listed in

Table 9. The AE has no fall-back implementation in the event that the Presentation Contexts are rejected by the SCP.

**Table 9: Proposed Presentation Contexts for AE Storage**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID		
Ophthalmic Photography 8-Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 SV1)	1.2.840.10008.1.2.4.70	SCU	None
Widefield Ophthalmic Photography Stereographic Projection Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.5	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 SV1)	1.2.840.10008.1.2.4.70	SCU	None
Ophthalmic Tomography Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.4	JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 SV1)	1.2.840.10008.1.2.4.70	SCU	None
Grayscale Softcopy Presentation State Storage SOP Class	1.2.840.10008.5.1.4.1.1.11.1	Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None
Stereometric Relationship Storage SOP Class	1.2.840.10008.5.1.4.1.1.77.1.5.3	Explicit VR Little Endian	1.2.840.10008.1.2.1	SCU	None

#### 4.4.1.3.1.3 SOP Specific Conformance for Storage SOP Classes

All Image & Presentation State Storage SOP Classes supported by the Storage AE exhibit the same behavior, except where stated, and are described together in this section.

The behavior of Storage AE when encountering status codes in a C-STORE response is summarized in the Table below:

**Table 10: Storage C-STORE Response Status Handling Behaviour**

Service Status	Further Meaning	Error Code	Behaviour
Success	Success	0000	The SCP has successfully stored the SOP Instance. If all SOP Instances in a send job have status success then the job is marked as complete.
*	*	Any other response	A failure is logged

A failed send job can be restarted by user interaction.

The contents of Ophthalmic Photography 8-bit Image Storage SOP Instances created by the OPTOS Suzuka software conform to the DICOM Ophthalmic Photography 8-bit Image IOD and are described in 7.1.1.1.

The contents of Widefield Ophthalmic Photography Stereographic Projection Image Storage SOP Instances created by the OPTOS Suzuka software conform to the DICOM Widefield Ophthalmic Photography Stereographic Projection Image IOD and are described in 7.1.1.2.

The contents of Ophthalmic Tomography Image Storage SOP Instances created by the OPTOS Suzuka software conform to the DICOM Ophthalmic Tomography Image IOD and are described in 7.1.1.3.

The contents of Grayscale Softcopy Presentation State Storage SOP Instances created by the OPTOS Suzuka software conform to the DICOM Grayscale Softcopy Presentation State IOD and are described in 7.1.1.4.

The contents of Stereometric Relationship Storage SOP Instances created by the OPTOS Suzuka software conform to the DICOM Stereometric Relationship IOD and are described in 7.1.1.4 .

#### **4.4.1.4 Association Acceptance Policy**

No Associations are accepted. This AE is SCU.

**4.4.1.4.1 Real-World Activity – for Modality Worklist**

Suzuka software initiates an association for DICOM C-FIND service.

**4.4.1.4.1.1 Associated Real-World Activity – for Modality Worklist**

The user performs a modality worklist query from the Optos device which is sent to the SCP via the Suzuka software. When Optos device queries for new worklist entries, Suzuka software opens an association to the Modality Worklist SCP. After the association has been established, a C-Find request with possible matching values is sent to the SCP. The provider then sends the matching worklist entries with C-Find responses to the Worklist SCU.

**4.4.1.4.1.2 Proposed Presentation Contexts – for Modality Worklist**

The presentation contexts that are proposed by Optos device and Suzuka software for the Modality Worklist C-Find operations are specified as follows:

**Table 11: Presentation Contexts for Modality Worklist**

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID		
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	Explicit VR Little Endian Implicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2 1.2.840.10008.1.2.2	SCU	None
Study Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.2.1	Explicit VR Little Endian Implicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2 1.2.840.10008.1.2.2	SCU	None
Patient Root Query/Retrieve Information Model - FIND	1.2.840.10008.5.1.4.1.2.1.1	Explicit VR Little Endian Implicit VR Little Endian Explicit VR Big Endian	1.2.840.10008.1.2.1 1.2.840.10008.1.2 1.2.840.10008.1.2.2	SCU	None

**4.4.1.4.1.3 SOP Specific Conformance for Modality Worklist**

The following Module descriptions contain the attributes, which are present in a C-FIND request messages sent by the Suzuka software to a remote AE. The descriptions also specify the attributes that are copied into the images.

Unexpected attributes returned in a C-FIND response are ignored.

The below tables should be read as follows:

Attribute Name: Attributes supported to build a Modality Worklist request.

Tag: DICOM tag for this attribute.

Required for MWL: These attributes are required for Optos devices and Suzuka software to be configured with Modality Worklist. { Yes=required; No=Not required }

Mapped into The Image: Whether this data is mapped into subsequently acquired images

□Notes: Clarification of this implementation’s use/treatment of this attribute

## SOP Common Module

**Table 12: SOP Common Module Attributes**

Attribute Name	Tag	Required for MWL	Mapped into the image	Notes
Specific Character Set	(0008,0005)	Yes	Yes	

## Patient List Query

**Table 13: Patient List Query Attributes**

Attribute Name	Tag	Required for MWL	Mapped into the image	Notes
Patient ID	(0010,0020)	No	Yes	Wild card matching is implemented for leading and trailing characters. Patient ID is displayed on the optos device screen
Patient Name	(0010,0010)	No	Yes	Matching is supported as follows: either no Patient's Name is supplied; the patient's last (family) name specified on the scanner's screen is supplied, or the Patient's first name for matching; this is user defined. Wild card matching is implemented for leading and trailing characters. Patient name is displayed on the optos device screen
Other Patient ID's	(0010,1000)	No	Yes	
Issuer of Patient ID	(0010,0021)	No		
Patient's Sex	(0010,0040)	Yes	Yes	Patient's Sex is displayed on the optos device screen
Patient's Birth Date	(0010,0030)	Yes	Yes	Patient's Date of birth is displayed on the optos device screen

## Study Query

All attributes as patient query, plus :

**Table 14: Study Query Attributes**

Attribute Name	Tag	Required for MWL	Mapped into the image	Notes
Study Instance UID	(0020,000D)	No		Unique identifier for the Study.
Modalities in Study	(0008,0061)	No		All of the distinct values used for Modality (0008,0060) in the Series of the Study.
Study ID	(0020,0010)	No		User or equipment generated

				Study identifier.
Accession Number	(0008,0050)	No		A RIS generated number that identifies the order of the study.
Study Date	(0008,0020)	No		Date the Study started.
Study Time	(0008,0030)	No		Time the Study started.
Study Description	(0008,1030)	No		Institution-generated description or classification of the Study (component) performed.
Number of Study Related Series	(0020,1206)	No		The number of series that match the Study level Query/Retrieve search criteria.
Number of Study Related Instances	(0020,1208)	No		The number of Composite Object Instances that match the Study level Query/Retrieve search criteria.

**Series Query**

**All attributes as patient query, plus :**

**Table 15: Series Query Attributes**

Attribute Name	Tag	Required for MWL	Mapped into the image	Notes
Study Instance UID	(0020,000D)	Yes		Unique identifier for the Study.
Series Instance UID	(0020,000E)	No		Unique identifier for the Series.
Series Number	(0020,0011)	No		A number that identifies the Series.
Series Description	(0008,103E)	No		Description of the Series.
Modality	(0008,0060)	No		Type of equipment that originally acquired the data used to create the images in this series.
Series Date	(0008,0021)	No		Date the Series started.
Series Time	(0008,0031)	No		Time the Series started.
Number of Series Related Instances	(0020,1209)	No		The number of Composite Object Instances in a Series that match the Series level Query/Retrieve search criteria.

**Work List Query**

**Table 16: Work List Query Attributes**

Attribute Name	Tag	Required for MWL	Mapped into the image	Notes
Patient ID	(0010,0020)	No	Yes	Wild card matching is implemented for leading and trailing characters. Patient ID is displayed on the optos device screen

Patient Name	(0010,0010)	No	Yes	Matching is supported as follows: either no Patient's Name is supplied; the patient's last (family) name specified on the scanner's screen is supplied, or the Patient's first name for matching; this is user defined. Wild card matching is implemented for leading and trailing characters. Patient name is displayed on the optos device screen
Other Patient ID's	(0010,1000)	No	Yes	
Issuer of Patient ID	(0010,0021)	No		
Patient's Sex	(0010,0040)	No	Yes	Patient's Sex is displayed on the optos device screen
Patient Weight	(0010,1030)	No		Weight of the Patient, in kilograms.
Patient's Birth Date	(0010,0030)	No	Yes	Patient's Date of birth is displayed on the optos device screen
Medical Alerts	(0010,2000)	No		Conditions to which medical staff should be alerted.
Pregnancy Status	(0010,21C0)	No		Describes pregnancy state of Patient.
Allergies	(0010,2110)	No		Description of prior reaction to contrast agents, or other patient allergies or adverse reactions.
Patient Comments	(0010,4000)	No		User-defined additional information about the Patient.
Special Needs	(0038,0050)	No		Medical and Social needs.
Patient State	(0038,0500)	No		Description of Patient state.
Current Patient Location	(0038,0300)	No		Describes the current known location of the Patient.
Institution Name	(0008,0080)	No		Institution where the equipment that produced the composite instances is located.
Admission ID	(0038,0010)	No		Visit and Admission are interchangeable here.
Accession Number	(0008,0050)	No		A RIS generated number that identifies the order of the study.
Referring Physician's Name	(0008,0090)	No	Yes	Referring Physician's Name is displayed on the optos device screen
Admitting Diagnoses	(0008,1080)	No		Description of the admitting

Description				diagnosis (diagnoses).
Requesting Physician	(0032,1032)	No		Name of the Physician who requested the Imaging Service Request.
Study Instance UID	(0020,000D)	Yes		Unique identifier for the Study.
Study Description	(0008,1030)	No		Institution-generated description or classification of the Study (component) performed.
Study ID	(0020,0010)	No		User or equipment generated Study identifier.
Reason for the Requested Procedure	(0040,1002)	No		Reason for requesting this procedure.
Study Date	(0008,0020)	No		Date the Study started.
Study Time	(0008,0030)	No		Time the Study started.
Requested Procedure ID	(0040,1001)	No		
Requested Procedure Description	(0032,1060)	No		
Requested Procedure Priority	(0040,1003)	No		Requested Procedure Type urgency.
Requested Procedure Code Sequence	(0032,1064)	No		
Referenced Study Sequence	(0008,1110)	No		
Procedure Code Sequence	(0008,1032)	No		A Sequence that conveys the type of procedure performed.

**Scheduled Procedure Step Module (for Worklist Query)**

**Table 17: Scheduled Procedure Step Module Attributes**

Attribute Name	Tag	Required for MWL	Mapped into the image	Notes
Scheduled Procedure Step Sequence	(0040,0100)	No		
> Scheduled Station AE Title	(0040,0001)	No		Suzuka software supports Broad query using scheduled AE Title
> Scheduled Station Name	(0040,0010)	No		Not used
> Scheduled Procedure Step Start Date	(0040,0002)	No		Suzuka software support Broad query using procedure start date
> Scheduled Procedure Step Start Time	(0040,0003)	No		Not used
> Modality	(0008,0060)	No	Yes	Matching is supported as follows: either no Modality is supplied (universal matching), or the scanner's Modality is supplied

				for matching; this is user defined.
> Scheduled Performing Physician's Name	(0040,0006)	No		Not used
> Scheduled Procedure Step Description	(0040,0007)	No		Not used
> Scheduled Protocol Code Sequence	(0040,0008)	No		Not used
> Scheduled Procedure Step Location	(0040,0011)	No		Not used
> Scheduled Procedure Step ID	(0040,0009)	No		Not used
> Requested Contrast Agent	(0032,1070)	No		Not used
> Pre-Medication	(0040,0012)	No		Not used
> Anatomical Orientation Type	(0010,2210)	No		Not used

## 4.5 Network interfaces

### 4.5.1 Physical Network Interface

The software operates in a Microsoft Windows environment, and as such is not associated with any specific physical networking hardware.

The software uses the Windows TCP/IP stack from the host machine.

### 4.5.2 Additional Protocols

The software does not make use of any additional protocols.

### 4.5.3 IPv4 and IPv6 Support

The software supports IPv4. The software does not support IPv6.

Where configuration permits and IP address to be specified, that IP address must be specified as an IPv4 address.

## 4.6 Configuration

### 4.6.1.1 Application Entity Title/Presentation Address Mapping

The local AE title is configurable through the configuration screen.

#### 4.6.1.1.1 Local Application Entity Titles

There is a configuration Screen for AE Titles, port numbers, and addresses. There are no configured defaults. The operator must configure values when enabling this feature.

**Table 18: AE Title Configuration Table**

Application Entity	Default AE Title	Default TCP/IP Port
Storage	No Default	No Default

### 4.6.1.2 Remote Application Entity Title / Presentation Address Mapping

The AE Title, host names and port numbers of remote applications are configured using the application configuration utility.

#### 4.6.1.2.1 Storage

The application configuration utility must be used to set the AE Titles, port-numbers, host-names and capabilities for the remote Storage SCPs. Associations will only be accepted from known AE Titles and associations from unknown AE Titles will be rejected (an AE Title is known if it can be selected within the Service/Installation Tool). Only a single remote Storage SCP can be defined.

### 4.6.1.3 Parameters

The following configuration parameters affect DICOM communication for the product:

**Table 19: Configuration Parameters**

Parameter	Configurable (Yes/No)	Default Value
<b>Storage Parameters</b>		
Support Widefield Ophthalmic Photography Stereographic Projection SOP Class	Yes	True
Character encoding for Image Storage SOP instances	Yes	ISO IR 192
Apply Lossless JPEG compression to images	Yes	True
Network name / IP address of the Store SCP (ImageProcessing_DicomServerHost)	Yes	localhost
TCP/IP port number for the Store SCP (ImageProcessing_DicomServerPort)	Yes	1100
Number of times to attempt a DIMSE C-STORE request. (ImageProcessingQueue_WorkerMaxNoOfHighRiskJob Assignments)	Yes	16 attempts

---

Calling AP title. (ImageProcessing_DicomClientAet)	Yes	NILAE
Called AP title. (ImageProcessing_DicomServerAet)	Yes	NILAE

## 5 Support for Extended Character Sets

### 5.1 Character sets in use

All Application Entities for the product support the following character sets:

- ISO IR 192 (UTF-8 encoding)
- ISO IR 100 (ISO 8859-1 : 1987 Latin alphabet No. 1 supplementary set)

#### 5.1.1 Character set configuration

The selection of character set for DICOM SOP instances is a user configurable setting within the software application. All DICOM SOP instances generated will use the same character set, as configured for the application.

By default, the application will use the ISO IR 192 character set.

## 6 Security

### 6.1 Conformance

The product does not claim conformance with the Secure Use Profiles listed in Part 15 of the DICOM specifications.

This product assumes that it operates within a secure environment.

## 7 Annexes

### 7.1 IOD Contents

#### 7.1.1 Created SOP Instances

Table 20 specifies the attributes of an Ophthalmic Photography 8-bit Image transmitted by OPTOS devices using the Suzuka storage application.

Table 21 specifies the attributes of a Widefield Ophthalmic Photography Stereographic Projection Image transmitted by OPTOS devices using the Suzuka storage application.

Table 22 specifies the attributes of an Ophthalmic Tomography Image transmitted by OPTOS devices using the Suzuka storage application.

Table 24 specifies the attributes of a Grayscale Softcopy Presentation State transmitted by OPTOS devices using the Suzuka storage application.

Abbreviations used for the Usage field in the tables are:

<b>VNAP</b>	Value Not Always Present (attribute sent zero length if no value is present)
<b>ANAP</b>	Attribute Not Always Present
<b>ALWAYS</b>	Always present with a value
<b>EMPTY</b>	Attribute is sent without a value

Abbreviations used for the source of the data values are:

<b>USER</b>	the attribute value source is from User input
<b>AUTO</b>	the attribute value is generated automatically
<b>CONFIG</b>	the attribute value source is a configurable parameter

Note

All date (DA) and time (TM) values are encoded in the local configured calendar and time. DateTime (DT) values will include an explicit timezone offset from UTC.

## 7.1.1.1 Ophthalmic Photography 8-Bit Image IOD

Table 20: IOD of Created Ophthalmic Photography 8-Bit Image Instances

IE	Module	Reference	Presence of Module
Patient	Patient	Table 26: Patient Module of Created SOP Instances	ALWAYS
Study	General Study	Table 27	ALWAYS
Series	General Series	Table 28	ALWAYS
	Ophthalmic Photography Series	Table 36	ALWAYS
Frame Of Reference	Frame Of Reference – <i>see Note 1</i>	Table 29	ALWAYS
	Synchronization	Table 30	ALWAYS
Equipment	General Equipment	Table 31	ALWAYS
Image	General Image	Table 37	ALWAYS
	Image Pixel	Table 38	ALWAYS
	Enhanced Contrast/Bolus	Table 33	Only present if contrast was administered
	Cine	N/A	Not used
	Multi-frame	Table 39	ALWAYS
	Ophthalmic Photography Image	Table 40	ALWAYS
	Ocular Region Imaged	Table 34	ALWAYS
	Ophthalmic Photography Acquisition Parameters	Table 41	ALWAYS
	Ophthalmic Photographic Parameters	Table 42	ALWAYS
	SOP Common	Table 35	ALWAYS
	Private Ophthalmic Photography Image – <i>See Note 1</i>	Table 43	ALWAYS

*Note 1* - Ophthalmic Photography 8-Bit Image SOP Instances conform to a Standard Extended SOP Class definition, with the addition of Frame of Reference and Private Ophthalmic Photography Image module attributes.

## 7.1.1.2 Widefield Ophthalmic Photography Stereographic Projection IOD

Table 21: IOD of Created Widefield Ophthalmic Photography SP Image Instances

IE	Module	Reference	Presence of Module
Patient	Patient	Table 26	ALWAYS
Study	General Study	Table 27	ALWAYS
Series	General Series	Table 28	ALWAYS
	Ophthalmic Photography Series	Table 44	ALWAYS

IE	Module	Reference	Presence of Module
Frame Of Reference	Frame Of Reference	Table 29	ALWAYS
	Synchronization	Table 30	ALWAYS
Equipment	General Equipment	Table 31	ALWAYS
	Enhanced General Equipment	Table 32	ALWAYS
Image	General Image	Table 45	ALWAYS
	Image Pixel	Table 46	ALWAYS
	Enhanced Contrast/Bolus	Table 33	Only present if contrast was administered
	Cine	N/A	Not used
	Multi-frame	Table 47	ALWAYS
	Ophthalmic Photography Image	Table 48	ALWAYS
	Widefield Ophthalmic Photography Stereographic Projection		ALWAYS
	Widefield Ophthalmic Photography Quality Rating	N/A	Not used
	Ocular Region Imaged	Table 34	ALWAYS
	Ophthalmic Photography Acquisition Parameters	Table 50	ALWAYS
	Ophthalmic Photographic Parameters	Table 51	ALWAYS
	ICC Profile	Table 52	Only if Photometric Interpretation (0028,0004) is not MONOCHROME2
	SOP Common	Table 35	ALWAYS
	Frame Extraction	N/A	Not used
	Private Ophthalmic Photography Image – see Note 2	Table 53	ALWAYS

Note 2 – Widefield Ophthalmic Photography Stereographic Projection Image SOP Instances conform to a Standard Extended SOP Class definition, with the addition of Private Ophthalmic Photography Image module attributes.

### 7.1.1.3 Ophthalmic Tomography Image IOD

Table 22: IOD of Created Ophthalmic Tomography Image Instances

IE	Module	Reference	Presence of Module
Patient	Patient	Table 26	ALWAYS
Study	General Study	Table 27	ALWAYS
Series	General Series	Table 28	ALWAYS
	Ophthalmic Tomography Series	Table 54	ALWAYS

IE	Module	Reference	Presence of Module
Frame Of Reference	Frame Of Reference	Table 29	ALWAYS
	Synchronization	Table 30	ALWAYS
Equipment	General Equipment	Table 31	ALWAYS
	Enhanced General Equipment	Table 32	ALWAYS
Image	General Image – see Note 3	Table 55	ALWAYS
	Image Pixel	Table 56	ALWAYS
	Enhanced Contrast/Bolus	N/A	Not used
	Multi-frame Functional Groups	Table 59	ALWAYS
	Multi-frame Dimension	Table 57	ALWAYS
	Acquisition Context	<b>Error! Reference source not found.</b>	ALWAYS
	Cardiac Synchronization	N/A	Not used
	Ophthalmic Tomography Image	Table 66	ALWAYS
	Ophthalmic Tomography Acquisition Parameters	Table 67	ALWAYS
	Ophthalmic Tomography Parameters	Table 68	ALWAYS
	Ocular Region Imaged	Table 34	ALWAYS
	SOP Common	Table 35	ALWAYS
	Private Ophthalmic Tomography Image – See Note 3	Table 69	ALWAYS

Note 3 - Ophthalmic Tomography Image SOP Instances conform to a Standard Extended SOP Class definition, with the addition of Private Ophthalmic Frame Location Functional Group Macro and Private Ophthalmic Tomography Image module attributes.

Table 23 below specifies the use of the Functional Group macros included in the Multi-Frame Functional Groups module for the Ophthalmic Tomography Image IOD.

**Table 23: Functional Group Macros Included in Ophthalmic Tomography Instances**

Functional Group Macro	Reference	Presence of Macro
Pixel Measures	Table 61	ALWAYS
Frame Content	Table 63	ALWAYS
Plane Position (Patient)	N/A	Not used
Plane Orientation (Patient)	N/A	Not used
Referenced Image	Table 60	ALWAYS
Derivation Image	N/A	Not used
Frame Anatomy	Table 62	ALWAYS
Cardiac Synchronization	N/A	Not used

---

Contrast/Bolus Usage	N/A	Not used
Ophthalmic Frame Location	Table 64	ALWAYS

## 7.1.1.4 Grayscale Softcopy Presentation State IOD

Table 24: IOD of Created Grayscale Softcopy Presentation State Instances

IE	Module	Reference	Presence of Module
Patient	Patient	Table 26	ALWAYS
Study	General Study	Table 27	ALWAYS
Series	General Series	Table 28	ALWAYS
	Presentation Series	Table 70	ALWAYS
Equipment	General Equipment	Table 31	ALWAYS
Presentation State	Presentation State Identification	Table 71	ALWAYS
	Presentation State Relationship	Table 72	ALWAYS
	Displayed Area	Table 73	ALWAYS
	Graphic Annotation	Table 74	Only if Graphic Annotations are present
	Graphic Layer		Only if Graphic Annotations are present
	Softcopy Presentation State LUT	Table 76	ALWAYS
	SOP Common	Table 35	ALWAYS

---

### 7.1.1.5 Stereometric Relationship IOD

**Table 25: IOD of Created Stereometric Relationship Storage Instances**

IE	Module	Reference	Presence of Module
Patient	Patient	Table 26	ALWAYS
Study	General Study	Table 27	ALWAYS
Series	General Series	Table 28	ALWAYS
	Stereometric Series	Table 77	ALWAYS
Equipment	General Equipment	Table 31	ALWAYS
Stereometric Relationship	Stereometric Relationship	Table 78	ALWAYS
	SOP Common	Table 35	ALWAYS

## 7.1.1.6 Common Modules

Table 26: Patient Module of Created SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Patient's Name	(0010,0010)	PN	Values supplied via User Input. Values will be a single group of 5 components (some possibly empty). Maximum 64 characters.	ALWAYS	USER
Patient ID	(0010,0020)	LO	Value supplied via User Input. Maximum 64 characters.	ALWAYS	USER
Patient's Birth Date	(0010,0030)	DA	From User Input.	ALWAYS	USER
Patient's Sex	(0010,0040)	CS	From User Input.	ALWAYS	USER
Patient Comments	(0010,4000)	LT	From User Input	ANAP	USER

Table 27: General Study Module of Created SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Study Instance UID	(0020,000D)	UI	Generated by device if not supplied by user.	ALWAYS	AUTO/ USER
Study Date	(0008,0020)	DA	Generated by device. <yyyymmdd>	ALWAYS	AUTO
Study Time	(0008,0030)	TM	Generated by device. <hhmmss.fff>	ALWAYS	AUTO
Referring Physician's Name	(0008,0090)	PN	Values supplied via User Input. Values will be a single group of 5 components (some possibly empty). Maximum 64 characters.	VNAP	USER
Study ID	(0020,0010)	SH	Zero Length	EMPTY	AUTO
Accession Number	(0008,0050)	SH	From User Input.	VNAP	USER
Study Description	(0008,1030)	LO	"Generic Eye Procedure 200 R"	ALWAYS	AUTO

Table 28: General Series Module of Created SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Modality	(0008,0060)	CS	"OP", "OPT" or "PR" depending upon the series type.	ALWAYS	AUTO
Series Instance UID	(0020,000E)	UI	Generated by the device.	ALWAYS	AUTO

Series Number	(0020,0011)	IS	Value generated by the device if Modality is OP or OPT. Value not present if Modality is PR.	VNAP	AUTO
Series Description	(0008,103E)	LO	Generated by the device, e.g. "RG OptomapPlus". Maximum 64 chars.	ALWAYS	AUTO
Operator's Name	(0008,1070)	PN	Name(s) of the operator(s) supporting the series	VNAP	AUTO
Body Part Examined	(0018,0015)	CS	If Modality is OP or OPT, value is "EYE". Value not present if Modality is PR.	VNAP	AUTO

**Table 29: Frame of Reference Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Frame Of Reference UID	(0020,0052)	UI	UID value generated by device. All images having a common Study Instance UID will share a common Frame of Reference UID.	ALWAYS	AUTO
Position Reference Indicator	(0020,1040)	LO	Generated from user input. "CORNEAL_VERTEX_L" or "CORNEAL_VERTEX_R"	ALWAYS	USER

**Table 30: Synchronisation Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Synchronization Trigger	(0018,106A)	CS	"NO TRIGGER"	ALWAYS	AUTO
Acquisition Time Synchronized	(0018,1800)	CS	"N"	ALWAYS	AUTO
Synchronization Frame Of Reference UID	(0020,0200)	UI	UID value generated by device. All images having a common Study Instance UID will share a common Synchronization Frame of Reference UID.	ALWAYS	AUTO

**Table 31: General Equipment Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Manufacturer	(0008,0070)	LO	"OPTOS"	ALWAYS	AUTO
Institution Name	(0008,0080)	LO	User input. Maximum 64 characters.	ALWAYS	USER
Station Name	(0008,1010)	SH	Generated by device, based on Device Serial Number. Maximum 16 characters.	ALWAYS	AUTO
Manufacturer's	(0008,1090)	LO	"P200T", "P200DTx", "P200TE" or	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
Model Name			"P200TxE"		
Device Serial Number	(0018,1000)	LO	Generated by device. Maximum 64 characters.	ALWAYS	AUTO
Software Version	(0018,1020)	LO	Current software versions of instrument. Maximum 64 characters.	ANAP	AUTO

**Table 32: Enhanced General Equipment Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Software Version	(0018,1020)	LO	Current software versions of instrument. Maximum 64 characters.	ALWAYS	AUTO

**Table 33: Enhanced Contrast/Bolus Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Contrast / Bolus Agent Sequence	(0018,0012)	SQ	Sequence contains one item:	ALWAYS	AUTO
> Code Value	(0008,0100)	SH	"C-B02CC" or "C-B0156"	ALWAYS	AUTO
> Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
> Code Meaning	(0008,0104)	LO	"Fluorescein" or "Indocyanine green"	ALWAYS	AUTO
> Contrast Bolus Agent Number	(0018,9337)	US	1	ALWAYS	AUTO
> Contrast / Bolus Administration Route Sequence	(0018,0014)	SQ	Sequence contains one item:	ALWAYS	AUTO
>> Code Value	(0008,0100)	SH	"G-D101"	ALWAYS	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>> Code Meaning	(0008,0104)	LO	"Intravenous route"	ALWAYS	AUTO
> Contrast/Bolus Ingredient Code Sequence	(0018,9338)	SQ	Zero Length	EMPTY	AUTO
> Contrast/Bolus Volume	(0018,1041)	DS	Zero Length	EMPTY	AUTO
> Contrast/Bolus Ingredient Concentration	(0018,1049)	DS	Zero Length	EMPTY	AUTO
> Contrast Administration Profile Sequence	(0018,9340)	SQ	Sequence contains one item:	ALWAYS	AUTO

>> Contrast Bolus Volume	(0018,1041)	DS	Zero Length	EMPTY	AUTO
>> Contrast Bolus Start Time	(0018,1042)	TM	Date and time generated by device. <yyyymmddhhmmss.fff&zzxx>	ALWAYS	AUTO

**Table 34: Ocular Region Imaged Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Anatomic Region Sequence	(0008,2218)	SQ	Sequence contains one item:	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"T-AA610"	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Retina"	ALWAYS	AUTO
Image Laterality	(0020,0062)	CS	From user input. "L" or "R".	ALWAYS	USER

**Table 35: SOP Common Module of Created SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Specific Character Set	(0008,0005)	CS	"ISO_IR 192" or "ISO_IR 100"	ALWAYS	CONFIG
SOP Class UID	(0008,0016)	UI	Generated by device. One of: 1.2.840.10008.5.1.4.1.1.77.1.5.1 1.2.840.10008.5.1.4.1.1.77.1.5.5 1.2.840.10008.5.1.4.1.1.77.1.5.4 1.2.840.10008.5.1.4.1.1.11.1	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Generated by device.	ALWAYS	AUTO

### 7.1.1.7 Ophthalmic Photography 8 Bit Image Modules

**Table 36: Ophthalmic Photography Series Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Modality	(0008,0060)	CS	"OP"	ALWAYS	AUTO

**Table 37: General Image Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Patient Orientation	(0020,0020)	CS	"L\F"	ALWAYS	AUTO
Image Comments	(0020,4000)	LT	User-defined comments about the image.	ANAP	USER

**Table 38: Image Pixel Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Rows	(0028,0010)	US	Rows in image.	ALWAYS	AUTO
Columns	(0028,0011)	US	Columns in image.	ALWAYS	AUTO
Bits Allocated	(0028,0100)	US	8	ALWAYS	AUTO
Bits Stored	(0028,0101)	US	8	ALWAYS	AUTO
High Bit	(0028,0102)	US	7	ALWAYS	AUTO
Pixel Data	(7FE0,0010)	OB	Image content. The Pixel data itself does not contain any burned-in annotation.	ALWAYS	AUTO

**Table 39: Multi-frame Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Number of Frames	(0028,0008)	US	1	ALWAYS	AUTO
Frame Increment Pointer	(0028,0029)	AT	(0008,002A) Acquisition DateTime	ALWAYS	AUTO

Table 40: Ophthalmic Photography Image Module of Created Ophth 8-bit SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Image Type	(0008,0008)	CS	Up to 8 fields separated by a backslash ('\') with defined values for each field as follows:- <b>Field #1</b> ORIGINAL <b>Field #2</b> PRIMARY <b>Field #3</b> <i>Empty field</i> <b>Field #4</b> <i>One of:</i> COLOR, RED, FA, ICG <b>Field #5 (category)</b> <i>One of:</i> OPTOMAP RG OPTOMAP AF OPTOMAPPLUS RG OPTOMAPPLUS AF OPTOMAPPLUS FA OPTOMAPPLUS ICG OCT INDEX IMAGE <b>Field #6</b> <i>Empty field</i> <b>Field #7</b> <i>Empty field</i> <b>Field #8</b> <i>If part of a stereometric pair, one of:</i> STEREO1 STEREO2	ALWAYS	AUTO
Instance Number	(0020,0013)	IS	Generated by device.	ALWAYS	AUTO
Samples Per Pixel	(0028,0002)	US	1 (monochrome) or 3 (composite).	ALWAYS	AUTO
Samples Per Pixel Used	(0028,0003)	US	Not present (monochrome) or 2 (composite).	ANAP	AUTO
Photometric Interpretation	(0028,0004)	CS	"MONOCHROME2" or "RGB".	ALWAYS	AUTO
Pixel Representation	(0028,0103)	US	0	ALWAYS	AUTO
Planar Configuration	(0028,0006)	US	Not present (monochrome) or 0 (composite).	ANAP	AUTO
Pixel Spacing	(0028,0030)	DS	Pixel spacing is mapped from device and image type. A pair of fixed-point numbers, to 4 decimal places: <f.ffff\f.ffff>	ALWAYS	AUTO
Content Time	(0008,0033)	TM	Time image pixel data creation started. <hhmmss.fff>	ALWAYS	AUTO
Content Date	(0008,0023)	DA	Date image pixel data creation started. <yyyymmdd>	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
Acquisition Date/Time	(0008,002A)	DT	Date and time generated by device. <yyyymmddhhmmss.fff&zzxx>	ALWAYS	AUTO
Lossy Image Compression	(0028,2110)	CS	00	ALWAYS	AUTO
Presentation LUT Shape	(2050,0020)	CS	Not used for Composite images. "IDENTITY" for Monochrome images.	ANAP	AUTO
Burned In Annotation	(0028,0301)	CS	"NO"	ALWAYS	AUTO

**Table 41: Ophthalmic Photography Acquisition Parameters Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Patient Eye Movement Commanded	(0022,0005)	CS	"NO" or "YES".	ALWAYS	USER
Patient Eye Movement Command Code Sequence	(0022,0006)	SQ	Sequence only present if Patient Eye Movement Commanded attribute has value of "YES".	ANAP	AUTO
>Code Value	(0008,0100)	SH	"R-404BF" or "R-404BC" or "R-404B6" or "R-404BD"	ALWAYS	USER
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Upward gaze" or "Left gaze" or "Downgaze" or "Right gaze"	ALWAYS	USER
Refractive State Sequence	(0022,001B)	SQ	Zero Length	EMPTY	AUTO
Emmetropic Magnification	(0022,000A)	FL	Zero Length	EMPTY	AUTO
Intra Ocular Pressure	(0022,000B)	FL	Zero Length	EMPTY	AUTO
Horizontal Field of View	(0022,000C)	FL	Zero length	EMPTY	AUTO
Pupil Dilated	(0022,000D)	CS	Zero Length	EMPTY	AUTO

**Table 42: Ophthalmic Photographic Parameters Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Acquisition Device Type Code Sequence	(0022,0015)	SQ	Sequence contains one item:	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"A-00E8A"	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Scanning Laser Ophthalmoscope"	ALWAYS	AUTO
Illumination Type Code Sequence	(0022,0016)	SQ	Zero Length	EMPTY	AUTO
Light Path Filter type stack code sequence	(0022,0017)	SQ	Zero Length	EMPTY	AUTO
Image Path Filter type stack code sequence	(0022,0018)	SQ	Zero Length	EMPTY	AUTO
Lenses Code Sequence	(0022,0019)	SQ	Zero Length	EMPTY	AUTO
Channel Description Code Sequence	(0022,001A)	SQ	Sequence of one or more of the following items:	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"G-A11A" or "G-A11E" or "G-A12F".	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Red", "Green" or "Blue"	ALWAYS	AUTO
Detector Type	(0018,7004)	CS	Zero Length	EMPTY	AUTO

**Table 43: Private Ophthalmic Photography Image Module of Created Ophth 8-bit SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Private Creator	(0021,0010)	LO	"OPTOS HEADER"	ALWAYS	AUTO
ProView	(0021,1000)	CS	"YES" or "NO"	ALWAYS	AUTO
SWL	(0021,1001)	CS	"NO"	ALWAYS	AUTO
Fovea Detection Algorithm Version	(0021,1010)	LO	Only present if ProView (0021,1000) is "YES". Maximum 64 chars.	ANAP	AUTO
Fovea Location	(0021,1011)	SL	Only present if ProView (0021,1000) is "YES". <x y> pixel co-ordinates; any negative value indicates an unknown location.	ANAP	AUTO
Fovea Confidence	(0021,1012)	FD	Only present if ProView (0021,1000) is "YES". Value between 0 and 1.	ANAP	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
Optic Disc Location	(0021,1013)	SL	Only present if ProView (0021,1000) is "YES". <x y> pixel co-ordinates; any negative value indicates an unknown location.	ANAP	AUTO
Optic Disc Confidence	(0021,1014)	FD	Only present if ProView (0021,1000) is "YES". Value 0 or greater.	ANAP	AUTO
Unprojected Fovea Location	(0021,1021)	SL	Only present if ProView (0021,1000) is "YES". Columns in unprojected image \ Rows in unprojected image; any negative value indicates an unknown location.	ANAP	AUTO
Unprojected Optic Disc Location	(0021,1022)	SL	Only present if ProView (0021,1000) is "YES". Columns in unprojected image \ Rows in unprojected image; any negative value indicates an unknown location.	ANAP	AUTO
Projection Algorithm Type	(0021,1023)	CS	Only present if ProView (0021,1000) is "YES". "FF" – fast fixed projection or "GP" – general projection	ANAP	AUTO
Projection Algorithm Table	(0021,1024)	LT	Only present if ProView (0021,1000) is "YES". Projection table name used in projection	ANAP	AUTO
Projection Algorithm Name	(0021,1025)	LT	Only present if ProView (0021,1000) is "YES". "OPTOS_FAST_FIXED" or "OPTOS_GENERAL_PROJECTION"	ANAP	AUTO
Projection Algorithm Version	(0021,1026)	LT	Only present if ProView (0021,1000) is "YES". Version of stereographic projection algorithm, e.g. "1.0.0"	ANAP	AUTO
Sharpening Applied	(0021,1031)	CS	Only present if Image Type category is one of: OPTOMAP RG, OPTOMAPPLUS RG Enumerated values: "YES" or "NO"	ANAP	AUTO
Sharpening Settings	(0021,1032)	DS	Only present if Image Type category is one of: OPTOMAP RG, OPTOMAPPLUS RG Repeated sequence of 4 fields per image channel (8 fields in total). Field values separated by a backslash ('\').	ANAP	AUTO

---

Attribute Name	Tag	VR	Value	Usage	Source
Optos Procedure Description	(0021, 1093)	LO	Value supplied via User Input. Maximum 64 characters.	VNAP	USER
Unprocessed Slo Image	(0021, 1098)	OB	Image binary. Only present if provided by the device.	ANAP	AUTO

### 7.1.1.8 Widefield Ophthalmic Photography Stereographic Projection Image Modules

**Table 44: Ophthalmic Photography Series Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Modality	(0008,0060)	CS	"OP"	ALWAYS	AUTO

**Table 45: General Image Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Patient Orientation	(0020,0020)	CS	"L\F"	ALWAYS	AUTO
Image Comments	(0020,4000)	LT	User-defined comments about the image.	ANAP	USER

**Table 46: Image Pixel Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Rows	(0028,0010)	US	Rows in image.	ALWAYS	AUTO
Columns	(0028,0011)	US	Columns in image.	ALWAYS	AUTO
Bits Allocated	(0028,0100)	US	8	ALWAYS	AUTO
Bits Stored	(0028,0101)	US	8	ALWAYS	AUTO
High Bit	(0028,0102)	US	7	ALWAYS	AUTO
Pixel Data	(7FE0,0010)	OB	Image content. The Pixel data itself does not contain any burned-in annotation.	ALWAYS	AUTO

**Table 47: Multi-frame Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Number of Frames	(0028,0008)	US	1	ALWAYS	AUTO
Frame Increment Pointer	(0028,0029)	AT	(0008,002A) Acquisition DateTime	ALWAYS	AUTO

**Table 48: Ophthalmic Photography Image Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Image Type	(0008,0008)	CS	Up to 8 fields separated by a backslash ('\') with defined values for each field as follows:- <b>Field #1</b> ORIGINAL <b>Field #2</b> PRIMARY <b>Field #3</b> <i>Empty field</i> <b>Field #4</b> <i>One of:</i> COLOR, RED, FA, ICG <b>Field #5 (category)</b> <i>One of:</i> OPTOMAP RG OPTOMAP AF OPTOMAPPLUS RG OPTOMAPPLUS AF OPTOMAPPLUS FA OPTOMAPPLUS ICG OCT INDEX IMAGE <b>Field #6</b> <i>Empty field</i> <b>Field #7</b> <i>Empty field</i> <b>Field #8</b> <i>If part of a stereometric pair, one of:</i> STEREO1 STEREO2	ALWAYS	AUTO
Instance Number	(0020,0013)	IS	Generated by device.	ALWAYS	AUTO
Samples Per Pixel	(0028,0002)	US	1 (monochrome) or 3 (composite).	ALWAYS	AUTO
Samples Per Pixel Used	(0028,0003)	US	Not present (monochrome) or 2 (composite).	ANAP	AUTO
Photometric Interpretation	(0028,0004)	CS	"MONOCHROME2" or "RGB".	ALWAYS	AUTO
Pixel Representation	(0028,0103)	US	0	ALWAYS	AUTO
Planar Configuration	(0028,0006)	US	Not present (monochrome) or 0 (composite).	ANAP	AUTO
Content Time	(0008,0033)	TM	Time image pixel data creation started. <hhmmss.fff>	ALWAYS	AUTO
Content Date	(0008,0023)	DA	Date image pixel data creation started. <yyyymmdd>	ALWAYS	AUTO
Acquisition Date/Time	(0008,002A)	DT	Date and time generated by device. <yyyymmddhhmmss.fff&zzxx>	ALWAYS	AUTO
Lossy Image Compression	(0028,2110)	CS	00	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
Presentation LUT Shape	(2050,0020)	CS	Not used for Composite images. "IDENTITY" for Monochrome images.	ANAP	AUTO
Burned In Annotation	(0028,0301)	CS	"NO"	ALWAYS	AUTO

**Table 49: Widefield Ophthalmic Photography Stereographic Projection Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Anatomic Region Sequence	(0008,2218)	SQ	Sequence containing a single item.	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"T-AA610"	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Retina"	ALWAYS	AUTO
> Anatomic Region Modifier Sequence	(0008,2218)	SQ	Sequence containing a single item.	ALWAYS	AUTO
>>Code Value	(0008,0100)	SH	"G-A100" or "G-A101"	ALWAYS	USER
>>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>>Code Meaning	(0008,0104)	LO	"Right" or "Left"	ALWAYS	USER
Transformation Algorithm Sequence	(0022,1513)	SQ	Sequence containing a single item:	ALWAYS	AUTO
> Algorithm Family Code Sequence	(0066,002F)	SQ	Sequence containing a single item:	ALWAYS	AUTO
>> Code Value	(0008,0100)	SH	"111791"	ALWAYS	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	"DCM"	ALWAYS	AUTO
>> Code Meaning	(0008,0104)	LO	"Spherical projection"	ALWAYS	AUTO
> Algorithm Name	(0066,0036)	LO	"OPTOS_FAST_FIXED" or "OPTOS_GENERAL_PROJECTION"	ALWAYS	AUTO
> Algorithm Version	(0066,0031)	LO	Version of stereographic projection algorithm, e.g. "1.0.0"	ALWAYS	AUTO
Ophthalmic Axial Length	(0022,1019)	FL	The axial length measurement used for the stereographic projection, in mm.	ALWAYS	AUTO
Ophthalmic Axial Length Method	(0022,1515)	CS	One of: "MEASURED", "ESTIMATED", "POPULATION"	ALWAYS	AUTO

X Coordinate Center Pixel View Angle	(0022,1528)	FL	The horizontal angle covered on the sphere by the center pixel of the projected image.	ALWAYS	AUTO
Y Coordinate Center Pixel View Angle	(0022,1529)	FL	The vertical angle covered on the sphere by the center pixel of the projected image.	ALWAYS	AUTO

**Table 50: Ophthalmic Photography Acquisition Parameters Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Patient Eye Movement Commanded	(0022,0005)	CS	"NO" or "YES".	ALWAYS	USER
Patient Eye Movement Command Code Sequence	(0022,0006)	SQ	Sequence only present if Patient Eye Movement Commanded attribute has value of "YES".	ANAP	AUTO
>Code Value	(0008,0100)	SH	"R-404BF" or "R-404BC" or "R-404B6" or "R-404BD"	ALWAYS	USER
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Upward gaze" or "Left gaze" or "Downgaze" or "Right gaze"	ALWAYS	USER
Refractive State Sequence	(0022,001B)	SQ	Zero Length	EMPTY	AUTO
Emmetropic Magnification	(0022,000A)	FL	Zero Length	EMPTY	AUTO
Intra Ocular Pressure	(0022,000B)	FL	Zero Length	EMPTY	AUTO
Horizontal Field of View	(0022,000C)	FL	Zero length	EMPTY	AUTO
Pupil Dilated	(0022,000D)	CS	Zero Length	EMPTY	AUTO

**Table 51: Ophthalmic Photographic Parameters Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Acquisition Device Type Code Sequence	(0022,0015)	SQ	N/A	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"A-00E8A"	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
>Code Meaning	(0008,0104)	LO	"Scanning Laser Ophthalmoscope"	ALWAYS	AUTO
Illumination Type Code Sequence	(0022,0016)	SQ	Zero Length	EMPTY	AUTO
Light Path Filter type stack code sequence	(0022,0017)	SQ	Zero Length	EMPTY	AUTO
Image Path Filter type stack code sequence	(0022,0018)	SQ	Zero Length	EMPTY	AUTO
Lenses Code Sequence	(0022,0019)	SQ	Zero Length	EMPTY	AUTO
Channel Description Code Sequence	(0022,001A)	SQ	Sequence of one or more of the following elements	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"G-A11A" or "G-A11E" or "G-A12F".	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Red", "Green" or "Blue"	ALWAYS	AUTO
Detector Type	(0018,7004)	CS	Zero Length	EMPTY	AUTO

Table 52: ICC Profile Module of Created Widefield Oph SP SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
ICC Profile	(0028,2000)	OB	An ICC Profile encoding the transformation of device-dependent color stored pixel values into PCS-Values	VNAP	AUTO

Table 53: Private Ophthalmic Photography Image Module of Created Widefield Oph SP SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Private Creator	(0021,0010)	LO	"OPTOS HEADER"	ALWAYS	AUTO
ProView	(0021,1000)	CS	"YES"	ALWAYS	AUTO
SWL	(0021,1001)	CS	"NO"	ALWAYS	AUTO
Fovea Detection Algorithm Version	(0021,1010)	LO	Only present if ProView (0021,1000) is "YES". Maximum 64 chars.	ANAP	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
Fovea Location	(0021,1011)	SL	Only present if ProView (0021,1000) is "YES". <x y> pixel co-ordinates; any negative value indicates an unknown location.	ANAP	AUTO
Fovea Confidence	(0021,1012)	FD	Only present if ProView (0021,1000) is "YES". Value between 0 and 1.	ANAP	AUTO
Optic Disc Location	(0021,1013)	SL	Only present if ProView (0021,1000) is "YES". <x y> pixel co-ordinates; any negative value indicates an unknown location.	ANAP	AUTO
Optic Disc Confidence	(0021,1014)	FD	Only present if ProView (0021,1000) is "YES". Value 0 or greater.	ANAP	AUTO
Unprojected Fovea Location	(0021,1021)	SL	Only present if ProView (0021,1000) is "YES". Columns in unprojected image \ Rows in unprojected image; any negative value indicates an unknown location.	ANAP	AUTO
Unprojected Optic Disc Location	(0021,1022)	SL	Only present if ProView (0021,1000) is "YES". Columns in unprojected image \ Rows in unprojected image; any negative value indicates an unknown location.	ANAP	AUTO
Projection Algorithm Type	(0021,1023)	CS	"FF" – fast fixed projection or "GP" – general projection	ANAP	AUTO
Projection Algorithm Table	(0021,1024)	LT	Projection table name used in projection	ANAP	AUTO
Projection Algorithm Name	(0021,1025)	LT	"OPTOS_FAST_FIXED" or "OPTOS_GENERAL_PROJECTION"	ANAP	AUTO
Projection Algorithm Version	(0021,1026)	LT	Version of stereographic projection algorithm, e.g. "1.0.0"	ANAP	AUTO
Sharpening Applied	(0021,1031)	CS	Only present if Image Type category is one of: OPTOMAP RG, OPTOMAPPLUS RG Enumerated values: "YES" or "NO"	ANAP	AUTO

---

Attribute Name	Tag	VR	Value	Usage	Source
Sharpening Settings	(0021,1032)	DS	Only present if Image Type category is one of: OPTOMAP RG, OPTOMAPPLUS RG Repeated sequence of 4 fields per image channel (8 fields in total). Field values separated by a backslash ('\').	ANAP	AUTO
Optos Procedure Description	(0021, 1093)	LO	Value supplied via User Input. Maximum 64 characters.	VNAP	USER
Unprocessed Slo Image	(0021, 1098)	OB	Image binary. Only present if provided by the device.	ANAP	AUTO

## 7.1.1.9 Ophthalmic Tomography Image Modules

**Table 54: Ophthalmic Tomography Series Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Modality	(0008,0060)	CS	"OPT"	ALWAYS	AUTO

**Table 55: General Image Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Patient Orientation	(0020,0020)	CS	"LVA"	ALWAYS	AUTO
Image Comments	(0020,4000)	LT	User-defined comments about the image.	ANAP	USER

**Table 56: Image Pixel Module of Created Widefield Oph SP SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Rows	(0028,0010)	US	Rows in each image frame.	ALWAYS	AUTO
Columns	(0028,0011)	US	Columns in each image frame.	ALWAYS	AUTO
Pixel Data	(7FE0,0010)	OB	Image content. The Pixel data itself does not contain any burned-in annotation.	ALWAYS	AUTO

**Table 57: Multi-frame Dimension Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Dimension Organization Sequence	(0020,9221)	SQ	Empty Group	EMPTY	AUTO
Dimension Index Sequence	(0020,9222)	SQ	Empty Group	EMPTY	AUTO

**Table 58: Acquisition Context Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Acquisition Context Sequence	(0040,0555)	SQ	Empty Group	EMPTY	AUTO

**Table 59: Multi-frame Functional Groups Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Shared Functional Groups Sequence	(5200,9229)	SQ	Sequence contains a single item:	EMPTY	AUTO
> Include Referenced Image Functional Group Macro			See Table 60	ALWAYS	N/A
> Include Frame Anatomy Functional Group Macro			See <b>Error! Reference source not found.</b>	ALWAYS	N/A
Per-frame Functional Groups Sequence	(5200,9230)	SQ	Sequence containing an item for every frame in the OCT image:	ALWAYS	AUTO
> Include Pixel Measures Functional Group Macro			See Table 61	ALWAYS	N/A
> Include Frame Content Functional Group Macro			See Table 63	ALWAYS	N/A
> Include Ophthalmic Frame Location Functional Group Macro			See Table 64	ALWAYS	N/A

**Table 60: Referenced Image Functional Group Macro of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Referenced Image Sequence	(0008,1140)	SQ	1 group	ALWAYS	AUTO
> Referenced SOP Class UID	(0008,1150)	UI	1.2.840.10008.5.1.4.1.1.77.1.5.1	ALWAYS	AUTO
> Referenced SOP Instance UID	(0008,1155)	UI	UID of SLO image instance	ALWAYS	AUTO
> Referenced Frame Number	(0008,1160)	IS	"1"	ALWAYS	AUTO
> Purpose of Reference Code Sequence	(0040,A170)	SQ	Sequence contains a single item:	ALWAYS	AUTO
>> Code Value	(0008,0100)	SH	"121311"	ALWAYS	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	"DCM"	ALWAYS	AUTO
>> Code Meaning	(0008,0104)	LO	"Localizer"	ALWAYS	AUTO
> Referenced SOP Class UID	(0008,1150)	UI	Only present if the image derived from or has created another image.  UID of referenced image class	ANAP	AUTO
> Referenced SOP Instance UID	(0008,1155)	UI	Only present if the image derived from or has created another image.  UID of referenced image instance	ANAP	AUTO

> Purpose of Reference Code Sequence	(0040,A170)	SQ	Only present if the image derived from or has created another image.  Sequence contains a single item:	ANAP	AUTO
>> Code Value	(0008,0100)	SH	Only present if the image derived from or has created another image.  "121338"	ANAP	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	Only present if the image derived from or has created another image.  "DCM"	ANAP	AUTO
>> Code Meaning	(0008,0104)	LO	Only present if the image derived from or has created another image.  "Anatomic image"	ANAP	AUTO
> Referenced SOP Class UID	(0008,1150)	UI	Only present if the image is a temporal successor to another image.  UID of referenced image class	ANAP	AUTO
> Referenced SOP Instance UID	(0008,1155)	UI	Only present if the image is a temporal successor to another image.  UID of referenced image instance	ANAP	AUTO
> Purpose of Reference Code Sequence	(0040,A170)	SQ	Only present if the image is a temporal successor to another image.  Sequence contains a single item:	ANAP	AUTO
>> Code Value	(0008,0100)	SH	Only present if the image is a temporal successor to another image.  "121348"	ANAP	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	Only present if the image is a temporal successor to another image.  "DCM"	ANAP	AUTO

>> Code Meaning	(0008,0104)	LO	Only present if the image derived from or has created another image.  "Temporal Predecessor"	ANAP	AUTO
-----------------	-------------	----	----------------------------------------------------------------------------------------------------	------	------

**Table 61: Pixel Measures Functional Group Macro of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Pixel Measures Sequence	(0028,9110)	SQ	Sequence contains a single item:	ALWAYS	AUTO
> Pixel Spacing	(0028,0030)	DS	Adjacent row spacing \ adjacent column spacing in mm. A pair of fixed point numbers to 4 decimal places: <f.fffff.ffff>	ALWAYS	AUTO

**Table 62: Frame Anatomy Functional Group Macro of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Frame Anatomy Sequence	(0020,9071)	SQ	Sequence contains a single item:	ALWAYS	AUTO
> Anatomic Region Sequence	(0008, 2218)	SQ	Sequence contains a single item:	ALWAYS	AUTO
>> Code Value	(0008,0100)	SH	"T-AA610"	ALWAYS	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	"SRT"	ALWAYS	AUTO
>> Code Meaning	(0008,0104)	LO	Retina	ALWAYS	AUTO
> Frame Laterality	(0020,9072)	CS	"R" = Right or OD, "L" =Left or OS	ALWAYS	USER

**Table 63: Frame Content Functional Group Macro of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Frame Content Sequence	(0020,9111)	SQ	Sequence contains a single item:	ALWAYS	AUTO
> Frame Acquisition Datetime	(0018,9074)	DT	<yyyymmddhhmmss.fff&zzxx>	ALWAYS	AUTO
> Frame Reference Datetime	(0018,9151)	DT	<yyyymmddhhmmss.fff&zzxx>	ALWAYS	AUTO
> Frame Acquisition Duration	(0018,9220)	FD	The actual amount of time [in milliseconds] that was used to acquire data for this frame	ALWAYS	AUTO

**Table 64: Ophthalmic Frame Location Functional Group Macro of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
----------------	-----	----	-------	-------	--------

Ophthalmic Frame Location Sequence	(0022,0031)	SQ	Sequence contains a single item:	ALWAYS	AUTO
> Referenced SOP Class UID	(0008,1150)	UI	1.2.840.10008.5.1.4.1.1.77.1.5.1	ALWAYS	AUTO
> Referenced SOP Instance UID	(0008,1155)	UI	UID of SLO image instance	ALWAYS	AUTO
> Referenced Frame Number	(0008,1160)	IS	"1"	ALWAYS	AUTO
> Purpose of Reference Code Sequence	(0040,A170)	SQ	Sequence contains a single item:	ALWAYS	AUTO
>>Code Value	(0008,0100)	SH	"121311"	ALWAYS	AUTO
>> Coding Scheme Designator	(0008,0102)	SH	"DCM"	ALWAYS	AUTO
>> Code Meaning	(0008,0104)	LO	"Localizer"	ALWAYS	AUTO
> Reference Coordinates	(0022,0032)	FL	Location on SLO Image. Row, Column values for the referenced image pixel coordinates.	ALWAYS	AUTO
> Ophthalmic Image Orientation	(0022,0039)	CS	"NONLINEAR" for RNFL otherwise "LINEAR"	ALWAYS	AUTO
> <i>Include Private Ophthalmic Frame Location Functional Group Macro</i>			<i>See table 55</i>	ALWAYS	N/A

**Table 65: Private Ophthalmic Frame Location Functional Group Macro of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Private Creator	(0021,0010)	LO	"OPTOS HEADER"	ALWAYS	AUTO
OCT Frame Location Confidence	(0021,1065)	IS	Percentage confidence in OCT frame location: Integer value between 0 and 100	ALWAYS	AUTO

**Table 66: Ophthalmic Tomography Image Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
----------------	-----	----	-------	-------	--------

Image Type	(0008,0008)	CS	Image Type String. 5 fields separated by a backslash ('\\') with defined values for each field as follows:- <b>Field #1</b> ORIGINAL <b>Field #2</b> PRIMARY <b>Field #3</b> POSTERIOR OCT or UWF OCT <b>Field #4</b> Empty field <b>Field #5 (category)</b> If field #3 is POSTERIOR OCT, then one of: ONH RNFL RETINA LINE RASTER If UWF OCT, then one of: NAV LINE NAV EXT LINE NAV VOLUME NAV HD VOLUME	ALWAYS	AUTO
Samples Per Pixel	(0028,0002)	US	1	ALWAYS	AUTO
Acquisition DateTime	(0008,002A)	DT	<yyyymmddhhss.fff&zzxx>	ALWAYS	AUTO
Acquisition Duration	(0018,9073)	FD	Varies by scan type	ALWAYS	AUTO
Acquisition Number	(0020,0012)	IS	"0"	ALWAYS	AUTO
Photometric Interpretation	(0028,0004)	CS	"MONOCHROME2"	ALWAYS	AUTO
Pixel Representation	(0028,0103)	US	0	ALWAYS	AUTO
Bits Allocated	(0028,0100)	US	8	ALWAYS	AUTO
Bits Stored	(0028,0101)	US	8	ALWAYS	AUTO
High Bit	(0028,0102)	US	7	ALWAYS	AUTO
Presentation LUT Shape	(2050,0020)	CS	"IDENTITY"	ALWAYS	AUTO
Lossy Image Compression	(0028,2110)	CS	"00"	ALWAYS	AUTO
Burned In Annotation	(0028,0301)	CS	"NO"	ALWAYS	AUTO
Concatenation Frame Offset Number	(0020,9228)	UL	0	ALWAYS	AUTO
In-concatenation Number	(0020,9162)	US	1	ALWAYS	AUTO
In-concatenation Total Number	(0020,9163)	US	1	ALWAYS	AUTO

Instance Number	(0020,0013)	IS	A number that identifies this instance	ALWAYS	AUTO
Content Date	(0008,0023)	DA	<yyyymmdd>	ALWAYS	AUTO
Content Time	(0008,0033)	TM	<hhmmss.fff>	ALWAYS	AUTO
Number of Frames	(0028,0008)	IS	Number of frames in a multi-frame image	ALWAYS	AUTO

**Table 67: Ophthalmic Tomography Acquisition Parameters Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Axial Length of the Eye	(0022,0030)	FL	Axial length of the eye in mm.	VNAP	AUTO
Refractive State Sequence	(0022,001B)	SQ	Empty group	EMPTY	AUTO
Emmetropic Magnification	(0022,000A)	FL	Zero length	EMPTY	AUTO
Intra Ocular Pressure	(0022,000B)	FL	Zero length	EMPTY	AUTO
Horizontal Field of View	(0022,000C)	FL	Zero length	EMPTY	AUTO
Pupil Dilated	(0022,000D)	CS	Zero length	EMPTY	AUTO

**Table 68: Ophthalmic Tomography Parameters Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Acquisition Device Type Code Sequence	(0022,0015)	SQ	Sequence contains a single item:	ALWAYS	AUTO
>Code Value	(0008,0100)	SH	"SRT"	ALWAYS	AUTO
>Coding Scheme Designator	(0008,0102)	SH	"A-00FBE"	ALWAYS	AUTO
>Code Meaning	(0008,0104)	LO	"Optical Coherence Tomography Scanner"	ALWAYS	AUTO
Light Path Filter Type Stack Code Sequence	(0022,0017)	SQ	Empty sequence	EMPTY	AUTO
Detector Type	(0018,7004)	CS	"CCD"	ALWAYS	AUTO
Illumination Wave Length	(0022,0055)	FL	Generated by device	ALWAYS	AUTO
Illumination Power	(0022,0056)	FL	Generated by device	ALWAYS	AUTO
Illumination Bandwidth	(0022,0057)	FL	Generated by device	ALWAYS	AUTO
Depth Spatial Resolution	(0022,0035)	FL	Generated by device	ALWAYS	AUTO
Maximum Depth Distortion	(0022,0036)	FL	Generated by device	ALWAYS	AUTO

Along-scan Spatial Resolution	(0022,0037)	FL	Generated by device	ALWAYS	AUTO
Maximum Along-scan Distortion	(0022,0038)	FL	Generated by device	ALWAYS	AUTO
Across-scan Spatial Resolution	(0022,0048)	FL	Generated by device	ALWAYS	AUTO
Maximum Across-scan Distortion	(0022,0049)	FL	Generated by device	ALWAYS	AUTO

**Table 69: Private Ophthalmic Tomography Image Module of Created OPT SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Private Creator	(0021,0010)	LO	"OPTOS HEADER"	ALWAYS	AUTO
RNFL Center Coordinates	(0021,1050)	FL	Only present if Image Type category is "RNFL". <x/y> pixel coordinates.	ANAP	AUTO
RNFL Start Coordinates	(0021,1051)	FL	Only present if Image Type category is "RNFL". <x/y> pixel coordinates.	ANAP	AUTO
RNFL Scan Direction	(0021,1052)	CS	Only present if Image Type category is "RNFL". One of: "CLOCKWISE" or "COUNTERCLOCKWISE"	ANAP	AUTO
OCT SNR Value	(0021,1060)	US	Integer value between 0 and 10	ALWAYS	AUTO
OCT Blink Detected Count	(0021,1061)	SH	"n/m" frames.	ANAP	AUTO
OCT Movement Detected Count	(0021,1062)	SH	"n/m" frames.	ANAP	AUTO
OCT Averages Count	(0021,1063)	SH	"n/m" frames.	ANAP	AUTO
OCT Scan Position Error	(0021,1064)	SH	Diagnostic message generated by device.	ANAP	AUTO
OCT Frame Location Confidence	(0021,1065)	IS	100	ANAP	AUTO
OCT Blink Detected	(0021,1066)	SH	Diagnostic message generated by device.	ANAP	AUTO
OCT Movement Detected	(0021,1067)	SH	Diagnostic message generated by device..	ANAP	AUTO
OCT Fovea Location	(0021,1071)	FL	Only present if segmentation analysis is enabled and the scan type is "Retina". <column\b-scan number> coordinates.	ANAP	AUTO

Attribute Name	Tag	VR	Value	Usage	Source
OCT Segmentation Outcome	(0021,1073)	SH	Only present if segmentation analysis is enabled. Diagnostic message generated by device.	ANAP	AUTO
OCT Segmentation Diagnostic	(0021,1074)	SH	Only present if the segmentation analysis is enabled but fails. Diagnostic message generated by device.	ANAP	AUTO
OCT Segmentation Disc Center Location	(0021, 107D)	FL	Only present if segmentation analysis is enabled and the scan type is "ONH". <column\b-scan number> coordinates.	ANAP	AUTO
OCT Segmentation Disc 3D Location	(0021, 107E)	FL	Only present if segmentation analysis is enabled and the scan type is "ONH". Multi-valued series of <column\row\b-scan number> coordinates indicating X/Y/Z of points describing the perimeter of the disc.	ANAP	AUTO
OCT Segmentation Cup 3D Location	(0021, 107F)	FL	Only present if segmentation analysis is enabled and the scan type is "ONH". Multi-valued series of <column\row\b-scan number> coordinates indicating X/Y/Z of points describing the perimeter of the cup.	ANAP	AUTO
Optos Procedure Description	(0021, 1093)	LO	Value supplied via User Input. Maximum 64 characters.	VNAP	USER
OCT SLO Reference Image	(0021, 1094)	OB	Image binary. Only present if provided by the device	ANAP	AUTO
OCT Scan Angle Degree	(0021, 1095)	FL	Only present if provided by the device Scan Rotation from reference image.	ANAP	AUTO
OCT Scan Center Location	(0021, 1096)	SL	Only present if provided by the device <x\y> pixel coordinates.	ANAP	AUTO
OCT SLO Thumbnail Image	(0021,1097)	OB	Image binary. Only present if provided by the device	ANAP	AUTO

---

Attribute Name	Tag	VR	Value	Usage	Source
OCT Segmentation Fovea Detection Confidence	(0021, 1099)	FL	Only present if segmentation analysis is enabled and the scan type is "Retina". Floating point number.	ANAP	AUTO
OCT Segmentation Fovea Detection Outcome	(0021, 109A)	SH	Only present if segmentation analysis is enabled and the scan type is "Retina". Diagnostic message generated by device.	ANAP	AUTO

## 7.1.1.10 Grayscale Softcopy Presentation State Modules

Table 70: Presentation Series Module of Created GSPS SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Modality	(0008,0060)	CS	"PR"	ALWAYS	AUTO
Series Description	(0008,103E)	LO	Derived from the Series Description of the series being presented, with a suffix of "(Presentation)". Maximum 64 chars.	ALWAYS	AUTO

Table 71: Presentation State Identification Module of Created GSPS SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Presentation Creation Date	(0070,0082)	DA	<yyyymmdd>	ALWAYS	AUTO
Presentation Creation Time	(0070,0083)	TM	<hhmmss.fff>	ALWAYS	AUTO
Instance Number	(0020,0013)	IS	One PR per study, so value is always "1"	ALWAYS	AUTO
Content Label	(0070,0080)	CS	"OPH SLO PR" or "RETINAL LAYERS"	ALWAYS	AUTO
Content Description	(0070,0081)	LO	"OPTOS RETINAL SLO IMAGE DISPLAY" or "OPTOS RETINAL LAYERS SEGMENTATION"	ALWAYS	AUTO
Content Creator's Name	(0070,0084)	PN	Zero length	EMPTY	AUTO

Table 72: Presentation State Relationship Module of Created GSPS SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Referenced Series Sequence	(0008,1115)	SQ	Sequence containing an item for every Series Instance UID, SOP Class UID combination that is presented by this PR.	ALWAYS	AUTO
> Series Instance UID	(0020,000E)	UI	UID of a series that is part of the same Study as this PR.	ALWAYS	AUTO
> Referenced Image Sequence	(0008,1140)	SQ	Sequence containing an item for every image within the referenced Series Instance UID that is presented this PR.	ALWAYS	AUTO
>> Referenced SOP Class UID	(0008,1150)	UI	SOP Class UID of presented image. All items within this sequence have the same SOP Class UID.	ALWAYS	AUTO
>> Referenced SOP Instance UID	(0008,1155)	UI	SOP Instance UID of presented image.	ALWAYS	AUTO

>> Referenced Frame Number	(0008,1160)	IS	Not present if all frames within the referenced SOP Instance UID are presented by this PR.	ANAP	AUTO
----------------------------	-------------	----	--------------------------------------------------------------------------------------------	------	------

**Table 73: Displayed Area Module of Created GSPS SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Displayed Area Selection Sequence	(0070,005A)	SQ	Sequence containing an item for every displayed area section. Every image and frame listed in the Presentation State Relationship Module is described:	ALWAYS	AUTO
> Referenced Image Sequence	(0008,1140)	SQ	Sequence containing an item for every image within some subset of the Presentation State Relationship Modul to which this displayed area selection applies:	ALWAYS	AUTO
>> Referenced SOP Class UID	(0008,1150)	UI	SOP Class UID of an image to which this displayed area selection applies.	ALWAYS	AUTO
>> Referenced SOP Instance UID	(0008,1155)	UI	SOP Instance UID of an image to which this displayed area selection applies.	ALWAYS	AUTO
>> Referenced Frame Number	(0008,1160)	IS	Not present if all frames within the referenced SOP Instance UID are included within this displayed area selection.	ANAP	AUTO
> Displayed Area Top Left Hand Corner	(0070,0052)	SL	Top-left pixel to be displayed. <x\y>	ALWAYS	AUTO
> Displayed Area Bottom Right Hand Corner	(0070,0053)	SL	Bottom-right pixel to be displayed. <x\y>	ALWAYS	AUTO
> Presentation Size Mode	(0070,0100)	CS	"SCALE TO FIT"	ALWAYS	AUTO

**Table 74: Graphic Annotation Module of Created GSPS SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Graphic Annotation Sequence	(0070,0001)	SQ	Sequence containing one or more items:	ALWAYS	AUTO
> Referenced Image Sequence	(0008,1140)	SQ	Sequence containing an item for every image within some subset of the Presentation State Relationship Modul to which this annotation selection applies:	ALWAYS	AUTO
>> Referenced SOP Class UID	(0008,1150)	UI	SOP Class UID of an image to which this annotation applies.	ALWAYS	AUTO

>> Referenced SOP Instance UID	(0008,1155)	UI	SOP Instance UID of an image to which this annotation applies.	ALWAYS	AUTO
>> Referenced Frame Number	(0008,1160)	IS	The frame to which this annotation applies. Not present if the annotation applies to all frames within the referenced image.	ANAP	AUTO
> Graphic Layer	(0070,0002)	CS	Label of the layer on which this annotation is to be displayed.	ALWAYS	AUTO
> Graphic Object Sequence	(0070,0009)	SQ	Sequence describing a graphic annotation. Contains one or more items.	ALWAYS	AUTO
>> Graphic Annotation Units	(0070,0005)	CS	"PIXEL"	ALWAYS	AUTO
>> Graphic Dimensions	(0070,0020)	US	2	ALWAYS	AUTO
>> Number of Graphic Points	(0070,0021)	US	Number of data points in this graphic.	ALWAYS	AUTO
>> Graphic Data	(0070,0022)	FL	Coordinates specifying the graphic annotation.	ALWAYS	AUTO
>> Graphic Type	(0070,0023)	CS	One of: "POINT", "POLYLINE"	ALWAYS	AUTO

**Table 75: Graphic Layer Module of Created GSPS SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Graphic Layer Sequence	(0070,0060)	SQ	Sequence containing one or more items:	ALWAYS	AUTO
> Graphic Layer	(0070,0002)	CS	Label used to identify the layer. Maximum 16 chars.	ALWAYS	AUTO
> Graphic Layer Order	(0070,0062)	IS	Integer indicating rendering order for layer.	ALWAYS	AUTO

**Table 76: Softcopy Presentation State LUT of Created GSPS SOP Instances**

Attribute Name	Tag	VR	Value	Usage	Source
Presentation LUT Shape	(2050,0020)	CS	"IDENTITY"	ALWAYS	AUTO

## 7.1.1.11 Stereometric Relationship Modules

Table 77: Stereometric Series Module of Created SRS SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Modality	(0008,0060)	CS	"SMR"	ALWAYS	AUTO

Table 78: Stereometric Relationship Module of Created SRS SOP Instances

Attribute Name	Tag	VR	Value	Usage	Source
Stereo Pairs Sequence	(0022,0020)	SQ	Sequence containing one or more items:	ALWAYS	AUTO
>Left Image Sequence	(0022,0021)	SQ	Sequence containing exactly one item	ALWAYS	AUTO
>> Referenced SOP Class UID	(0008,1150)	UI	Uniquely identifies the referenced SOP Class	ALWAYS	AUTO
>> Referenced SOP Instance UID	(0008,1155)	UI	Uniquely identifies the referenced SOP Instance	ALWAYS	AUTO
>Right Image Sequence	(0022,0022)	SQ	Sequence containing exactly one item	ALWAYS	AUTO
>> Referenced SOP Class UID	(0008,1150)	UI	Uniquely identifies the referenced SOP Class	ALWAYS	AUTO
>> Referenced SOP Instance UID	(0008,1155)	UI	Uniquely identifies the referenced SOP Instance	ALWAYS	AUTO

## 7.2 Usage of Attributes from received IOD's

Not applicable.

## 7.3 Attribute Mapping

Not applicable.

## 7.4 Coerced / Modified Fields

Not applicable.

## 7.5 Data Dictionary of Private Attributes

The Private Attributes added to created SOP Instances are listed in the table below. OPTOS Suzuka software reserves blocks of private attributes in group 0021.

**Table 79: Private Data Elements of Created SOP Instances**

Attribute Name	Tag	VR	VM	Attribute Description
Creator	(0021,0010)	LO	1	"OPTOS HEADER"
Preview	(0021,1000)	CS	1	"YES"/"NO"
SWL	(0021,1001)	CS	1	"YES"/"NO"
Fovea Detection Algorithm Version	(0021,1010)	LO	1	xxx.xx.xx
Fovea Location	(0021,1011)	SL	2	X and Y pixel coordinates
Fovea Confidence	(0021,1012)	FD	1	Between 0 and 1
Optic Disc Location	(0021,1013)	SL	2	X and Y pixel coordinates
Optic Disc Confidence	(0021,1014)	FD	1	Greater than 0
Unprojected Fovea Location	(0021,1021)	SL	2	X and Y pixel coordinates
Unprojected Optic Disc Location	(0021,1022)	SL	2	X and Y pixel coordinates
Projection Algorithm Type	(0021,1023)	CS	1	"FF"/"GP"
Projection Algorithm Table	(0021,1024)	LT	1	
Projection Algorithm Name	(0021,1025)	LT	1	
Projection Algorithm Version	(0021,1026)	LT	1	xxx.xx.xx
Sharpening Applied	(0021,1031)	CS	1	"YES"/"NO"
Sharpening Settings	(0021,1032)	DS	4-4n	Image sharpening settings. 4 values per image channel.
Auto Contrast	(0021,1040)	CS	1	"MANUAL"/"AUTO"
Detector Gain	(0021,1041)	IS	1	Integer value greater than 0

Detector Contrast	(0021,1042)	IS	1	Integer value greater than 0
RNFL Center Coordinates	(0021,1050)	FL	2	X and Y pixel coordinates
RNFL Start Coordinates	(0021,1051)	FL	2	X and Y pixel coordinates
RNFL Scan Direction	(0021,1052)	CS	1	"CLOCKWISE" or "COUNTERCLOCKWISE"
OCT SNR Value	(0021,1060)	US	1	Integer value between 0 and 10
OCT Blink Detected Count	(0021,1061)	SH	1	"n/m" where n is the number of frames in which blink was detected and m is the maximum number of frames in the scan. This field is populated if modality is OPT. In case, where n is 0, an empty value shall be encoded.
OCT Movement Detected Count	(0021,1062)	SH	1	"n/m" where n is the number of frames in which movement was detected and m is the maximum number of frames in the scan. This field is populated if modality is OPT. In case, where n is 0, an empty value shall be encoded.
OCT Averages Count	(0021,1063)	SH	1	"n/m" where n is the number of averaged frames and m is the maximum number of frames in the scan. This field is populated if modality is OPT.
OCT Scan Position Error	(0021,1064)	SH	1	"ERROR"/""
OCT Frame Location Confidence	(0021,1065)	IS	1	Percentage confidence in OCT frame location: Integer value between 0 and 100
OCT Blink Detected	(0021,1066)	SH	1	"WARNING"/""
OCT Movement Detected	(0021,1067)	SH	1	"WARNING"/""
OCT Fovea Location	(0021,1071)	FL	2	X pixel coordinate and the OCT frame number. This field is populated if the scan type is Retina.
OCT Segmentation Outcome	(0021,1073)	SH	1	"Seg Success"/"" or "Seg Fail"/""
OCT Segmentation Diagnostic	(0021,1074)	SH	1	"n" where n is the diagnostic code. This field is populated if the segmentation analysis fails.
OCT Segmentation Disc Center Location	(0021,107D)	FL	2	Only present if segmentation analysis is enabled and the scan type is "ONH". <column\b-scan number> coordinates.
OCT Segmentation Disc 3D Location	(0021,107E)	FL	3-3n	Multi-valued series of <column\row\b-scan number> coordinates indicating X/Y/Z of points describing the perimeter of the disc.

OCT Segmentation Cup 3D Location	(0021,107F)	FL	3-3n	Multi-valued series of <column\row\b-scan number> coordinates indicating X/Y/Z of points describing the perimeter of the cup.
Optos Procedure Description	(0021,1093)	LO	1	Optional free form text to label a study e.g. "Age Related Eye Disease Study".
OCT Slo Reference Image	(0021,1094)	OB	1	Image to be supplied when details of the OCT are requested by the scanhead
OCT Scan Angle Degree	(0021,1095)	FL	1	Angle to be supplied when details of the OCT are requested by the scanhead
OCT Scan Center Location Pixels	(0021,1096)	SL	1	Location to be supplied when details of the OCT are requested by the scanhead
Thumbnail Image	(0021,1097)	OB	1	Image to be supplied when details of the OCT are requested by the scanhead
Unprojected Slo Image	(0021,1098)	OB	1	Unprojected Slo image that was used for OCT planning.
OCT Segmentation Retina Fovea Detection Confidence	(0021, 1099)	FL	1	Confidence in detection of the OCT fovea location.
OCT Segmentation Retina Fovea Detection Outcome	(0021,109A)	SH	1	"Fovea Found"/"" or "Fovea Not Found"/""

## 7.6 Coded terminology and Templates

Not Applicable.

## 7.7 Grayscale Image Consistency

Not applicable.

## 7.8 Standard/Extended/Specialized/Private SOP Classes

The product provides Standard Extended Conformance to all supported SOP Classes, through the inclusion of additional Type 3 Standard Elements and Private Data Elements.

The extensions are defined in Section 7.1.1.

## 7.9 Private Transfer Syntaxes

Not applicable